

# Effects of a novel brace for treatment of the patients with adolescent idiopathic scoliosis; Good compliance leads to satisfactory results.

Kai HIRATA<sup>1)</sup> Osamu SHIRADO<sup>2)</sup> Masanori ABE<sup>3)</sup> Takuya KIGAWA<sup>1)</sup> Tatsuya ENDO<sup>2)</sup> Masumi IWABUCHI<sup>2)</sup>

1) Nippon Sigmax Co., Ltd

2) Aizu Medical Center Orthopaedics & Spinal Surgery, Fukushima Medical University

3) Tohoku Prosthetic and Orthotic Co., Ltd



## Purpose

Many studies have shown that brace treatment is the only effective way to treat adolescent idiopathic scoliosis (AIS) conservatively. However, the compliance for brace treatment is low. We developed a novel brace and monitoring system to solve this problem.

**The purpose of this study was to investigate the middle-term effects of a novel brace from both sides of the correction rate and wearing time.**

## Introduction of AMEC brace and monitoring system

### 【AMEC Brace】

Based on the principle of 3 point-fixation.



- ① Trunk Cuff : Made of FRTP
- ② Pad : Dial adjustment system


Concept : Light, Low profile, Easy to wear

### 【Monitoring system】

It consists of a temperature sensor, a pressure sensor, and a data logger.

Measurement interval : every **10** minutes

Both data can be saved for **300** days.



Brace compliance :

→『no wear』, 『wear』, 『effective wear』

## Methods

The subjects are **14** patients with AIS.  
The AMEC brace treatment was started with the consent of patients and guardians.  
The follow-up period is 1 year.  
We observed the brace compliance using the monitoring system for 1 month.

### 【Inclusion criteria】

- 10~15 years
- Cobb angle of 25~40 degrees
- Risser sign less than 3
- Within 2 years after menarche


## Results and Discussions

There were **8** lumbar and **6** thoracic curves.

	Initial correction rate	Correction rate after 1 year	wearing time after 1 year
Lumbar (N=8)	43.7 %	<b>42.9 %</b>	<b>16.3</b> hours/day
Thoracic (N=6)	29.0 %	22.6 %	9.3 hours/day

The AMEC brace was enough effective to treat the patients with AIS with high compliance. It was especially effective for the lumbar curve. Thoracic type modified brace design.


### Correction rate (Typical)



#### 【Lumbar curve】

- 14years old
- Female
- Risser sign 3

Correction rate  
→ **41.2%**




#### 【Thoracic curve】


- 11years old
- Female
- Risser sign 1

Correction rate  
→ **33.3%**

### Wearing time (Typical)



Wearing time: **21.2** hours/day



Effective wearing time: **13.9** hours/day

As in the past study(Weinstein, 2013), correction rate after 1 year was high for lumbar curve with long wearing time.

## Conclusion

**We conclude that good compliance for brace treatment leads to satisfactory results.**